

Product Datasheet

Lactoferrin Antibody HM810

Product Name	Lactoferrin Antibody
Catalog Number	HM810
Species/Host	Mouse
Tested Applications	LF, ELISA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with bovine lactoferrin, human/pig/bovine hemoglobin, human calprotectin, human/pig/bovine transferrin
Matched Pairs	HM811 (capture)
Preservatives	0.09% Sodium Azide