

Product Datasheet

D-dimer Antibody HM812

Product Name	D-dimer Antibody
Catalog Number	HM812
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF, TIA
Immunogen	Native D-dimer
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4 ± 0.2
Matched Pairs	HM812 (self-pairs)
Preservatives	0.05% Proclin 300