

## Product Datasheet

### Hemoglobin A1c (HbA1c) Antibody HM813

<b>Product Name</b>	Hemoglobin A1c (HbA1c) Antibody
<b>Catalog Number</b>	HM813
<b>Description</b>	Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	WB, ELISA, LF, TIA
<b>Immunogen</b>	Native HbA1c with peptides
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 90 %
<b>Buffer</b>	Phosphate Buffered Saline, pH 7.4 ± 0.2
<b>Matched Pairs</b>	HM814 (capture)
<b>Preservatives</b>	0.05% Proclin 300