

Product Datasheet

Hemoglobin A1c (HbA1c) Antibody HM814

Product Name	Hemoglobin A1c (HbA1c) Antibody
Catalog Number	HM814
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	WB, LF, ELISA
Immunogen	Native HbA1c with peptides
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM813 (detection)
Preservatives	0.05% Proclin 300