

Product Datasheet

Brain Natriuretic Peptide (BNP) Antibody HM818

Product Name	Brain Natriuretic Peptide (BNP) Antibody
Catalog Number	HM818
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, CLIA
Immunogen	Recombinant BNP
Alternative Names	B-type Natriuretic Peptide
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2a
Clonality	Monoclonal
Epitope	26-32 a.a.
Purity	> 90 %
Buffer	Phosphate Buffered Saline, pH 7.4 ± 0.2
Matched Pairs	HM819 (capture)
Preservatives	0.05% Proclin 300