

Product Datasheet

Malaria (pf-HRP-2) Antibody HM824

Product Name	Malaria (pf-HRP-2) Antibody
Catalog Number	HM824
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant Plasmodium falciparum HRP-2
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4 with 150 mM NaCl
Specificity	Plasmodium falciparum (malaria) HRP-2
Matched Pairs	HM825 (detection)
Preservatives	0.1% Proclin 300