

Product Datasheet

Human IL-2 ELISA Set 851.500

| | |
|---------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Name | Human IL-2 ELISA Set |
| Catalog Number | 851.500 |
| Description | Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition. |
| Range | 31.2-1000 pg/ml |
| Sensitivity | 7 pg/ml |
| Tested Applications | ELISA |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Cross Reactivity | No cross reactivity with other human cytokines |
| Specificity | Recognizes both natural and recombinant human IL-2 |
| Sample Types | Serum, Plasma, Cell culture supernatant |