

Product Datasheet

Malaria (pf-HRP-2) Antibody HM825

Product Name	Malaria (pf-HRP-2) Antibody
Catalog Number	HM825
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Recombinant Plasmodium falciparum HRP-2
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	-20°C. Avoid repeated freezing and thawing.
Clonality	Monoclonal
Purity	≥90% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4 with 150 mM NaCl
Specificity	Plasmodium falciparum (malaria) HRP-2
Matched Pairs	HM824 (capture)
Preservatives	0.1% Proclin 300