

Product Datasheet

Malaria (pf-HRP-2) Antibody HM825

Product Name	Malaria (pf-HRP-2) Antibody
Catalog Number	HM825
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant Plasmodium falciparum HRP-2
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	-20 °C. Avoid repeated freezing and thawing.
Clonality	Monoclonal
Purity	≥ 90 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4 with 150 mM NaCl
Specificity	Plasmodium falciparum (malaria) HRP-2
Matched Pairs	HM824 (capture)
Preservatives	0.1% Proclin 300