

Product Datasheet

Foot and Mouth Disease Virus (FMDV-Type O) Antibody HM833

Product Name	Foot and Mouth Disease Virus (FMDV-Type O) Antibody
Catalog Number	HM833
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Native FMDV
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	>90% by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Specificity	FMDV 3AB nonstructural protein
Matched Pairs	HM834, HM835, HM836 (detection)
Preservatives	0.05% Sodium Azide