

Product Datasheet

Haptoglobin, Feline Antibody HM837

Product Name	Haptoglobin, Feline Antibody
Catalog Number	HM837
Species/Host	Mouse
Tested Applications	LF, ELISA, WB
Immunogen	Feline Haptoglobin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	2-8°C
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM838 or HM839 (detection)
Preservatives	0.1% Sodium Azide