

## Product Datasheet

### Haptoglobin, Feline Antibody HM838

<b>Product Name</b>	Haptoglobin, Feline Antibody
<b>Catalog Number</b>	HM838
<b>Description</b>	Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	WB, ELISA, LF
<b>Immunogen</b>	Feline Haptoglobin
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	+2-8 °C
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 90 %
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.2
<b>Matched Pairs</b>	HM837 (capture)
<b>Preservatives</b>	0.05% Proclin 300