

Product Datasheet

Serum Amyloid A (SAA), Horse Antibody HM840

Product Name	Serum Amyloid A (SAA), Horse Antibody
Catalog Number	HM840
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	WB, ELISA, LF
Immunogen	Horse SAA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (< 2 weeks): +2-8 °C. Long Term: -20 °C
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM840 (self-pairs)
Preservatives	0.1% Sodium Azide