

Product Datasheet

Bovine Growth Hormone (rBGH) Antibody HM846

Product Name	Bovine Growth Hormone (rBGH) Antibody
Catalog Number	HM846
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	rBGH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	Pituitary Bovine Growth Hormone (pBGH), Native Bovine Growth Hormone (nBGH)
Specificity	Recombinant Bovine Growth Hormone
Matched Pairs	HM845 or HM847 (detection)
Preservatives	0.05% Sodium Azide