

Product Datasheet

Creatine Kinase MB (CK-MB) Antibody HM848

Product Name	Creatine Kinase MB (CK-MB) Antibody
Catalog Number	HM848
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF, TIA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Matched Pairs	HM849 (detection)
Preservatives	0.1% Proclin 300