

Product Datasheet

Creatine Kinase MB (CK-MB) Antibody HM849

| | |
|----------------------------|--|
| Product Name | Creatine Kinase MB (CK-MB) Antibody |
| Catalog Number | HM849 |
| Description | Purified Monoclonal Antibody |
| Species/Host | Mouse |
| Tested Applications | ELISA, LF, TIA |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | > 90 % |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.4 |
| Matched Pairs | HM848 (capture) |
| Preservatives | 0.1% Proclin 300 |