

Product Datasheet

Calicivirus, Feline Antibody HM854

Product Name	Calicivirus, Feline Antibody
Catalog Number	HM854
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Native Feline Calicivirus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	Phosphate Buffered Saline, pH 7.0
Matched Pairs	HM855 (capture)
Preservatives	0.09% Sodium Azide