

Product Datasheet

C-Reactive Protein, Canine (cCRP) Antibody HM857

Product Name	C-Reactive Protein, Canine (cCRP) Antibody
Catalog Number	HM857
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Native Canine CRP
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG
Clonality	Monoclonal
Purity	>95% by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.0
Matched Pairs	HM862 (detection)
Preservatives	0.09% Sodium Azide