

Product Datasheet

Parvovirus, Canine Antibody HM860

Product Name	Parvovirus, Canine Antibody
Catalog Number	HM860
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	WB, ELISA, LF
Immunogen	Canine Parvovirus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≥ 2 weeks) +2-8 °C. Long Term: -20 °C or below
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Matched Pairs	HM861 (detection)
Preservatives	0.09% Sodium Azide