

Product Datasheet

Pepsinogen II (PGII) Antibody HM876

Product Name	Pepsinogen II (PGII) Antibody
Catalog Number	HM876
Description	Purified Monoclonal Antibody
Species/Host	CHO Cells
Tested Applications	ELISA, CLIA, FIA, TIA
Immunogen	Recombinant pepsinogen II
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	PGII
Matched Pairs	HM877 (detection)
Preservatives	0.1% Proclin 300