

## Product Datasheet

### Pepsinogen II (PGII) Antibody HM877

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Pepsinogen II (PGII) Antibody  |
| <b>Catalog Number</b>      | HM877  |
| <b>Description</b>         | Purified Monoclonal Antibody   |
| <b>Species/Host</b>        | CHO Cells  |
| <b>Tested Applications</b> | ELISA, CLIA, FIA, TIA  |
| <b>Immunogen</b>           | Recombinant pepsinogen II  |
| <b>Brand</b>               | EastCoast Bio  |
| <b>Form/Appearance</b>     | Purified Monoclonal Antibody   |
| <b>Storage</b>             | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| <b>Isotype</b>             | IgG1   |
| <b>Clonality</b>           | Monoclonal   |
| <b>Purity</b>              | > 90 %   |
| <b>Buffer</b>              | 10 mM Phosphate Buffered Saline, pH 6.0                                      |
| <b>Specificity</b>         | PGII   |
| <b>Matched Pairs</b>       | HM876 (capture)  |
| <b>Preservatives</b>       | 0.1% Proclin 300   |