

# Product Datasheet

## Human IL-4 ELISA Set 851.510.020

Product Name	Human IL-4 ELISA Set
Catalog Number	851.510.020
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	1.1-35 pg/ml
Sensitivity	0.7 pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes both natural and recombinant human IL-4
Sample Types	Serum, Plasma, Cell culture supernatant