

# Product Datasheet

## Human IL-6 ELISA Set 851.520.005

|                     |  |
|---------------------|--|
| Product Name        | Human IL-6 ELISA Set   |
| Catalog Number      | 851.520.005  |
| Description         | Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition. |
| Range               | 6.25-200 pg/ml   |
| Sensitivity         | 2 pg/ml  |
| Tested Applications | ELISA  |
| Alternative Names   | NULL   |
| Brand               | Diaclone   |
| Storage             | Store at +4°C  |
| Cross Reactivity    | No cross reactivity with other human cytokines. No interference with soluble receptor gp80.  |
| Specificity         | Recognizes both natural and recombinant human IL-6   |
| Sample Types        | Serum, Plasma, Cell culture supernatant  |