

Product Datasheet

Human IL-6 ELISA Set 851.520.010

Product Name	Human IL-6 ELISA Set
Catalog Number	851.520.010
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	6.25-200 pg/ml
Sensitivity	2 pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. No interference with soluble receptor gp80.
Specificity	Recognizes both natural and recombinant human IL-6
Sample Types	Serum, Plasma, Cell culture supernatant