

Product Datasheet

Treponema pallidum (T. pallidum) Antibody HM894

Product Name	Treponema pallidum (T. pallidum) Antibody
Catalog Number	HM894
Species/Host	Mouse
Tested Applications	ELISA, WB, IFA
Alternative Names	Syphilis
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline pH 7.2
Cross Reactivity	Treponema reiter
Specificity	T. pallidum ~32 kDa protein
Preservatives	0.1% Sodium Azide