

# Product Datasheet

## Troponin I (41-49) Antibody HM897

|                     |                                                                              |
|---------------------|------------------------------------------------------------------------------|
| Product Name        | Troponin I (41-49) Antibody                                                  |
| Catalog Number      | HM897                                                                        |
| Description         | Purified Monoclonal Antibody                                                 |
| Species/Host        | CHO Cells                                                                    |
| Tested Applications | LF                                                                           |
| Brand               | EastCoast Bio                                                                |
| Form/Appearance     | Purified Monoclonal Antibody                                                 |
| Storage             | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Clonality           | Monoclonal                                                                   |
| Epitope             | 41-49 a.a.                                                                   |
| Purity              | > 90 %                                                                       |
| Buffer              | 10 mM Phosphate Buffered Saline, pH 7.4                                      |
| Specificity         | Troponin I                                                                   |
| Matched Pairs       | HM898 (capture) and HM896 (detection)                                        |
| Preservatives       | 0.1% Proclin 300                                                             |