

## Product Datasheet

### NT-proBNP (43-50) Antibody HM899

<b>Product Name</b>	NT-proBNP (43-50) Antibody
<b>Catalog Number</b>	HM899
<b>Description</b>	Purified Monoclonal Antibody
<b>Species/Host</b>	Sheep
<b>Tested Applications</b>	CLIA, LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	Short Term ( $\leq 2$ weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 90 %
<b>Buffer</b>	10 mM Phosphate Buffer, pH 7.4 $\pm$ 0.2 with 0.15M NaCl
<b>Specificity</b>	NT-proBNP (a.a. 43-50)
<b>Matched Pairs</b>	HM148 or HM210 (detection)
<b>Preservatives</b>	0.1% Proclin 300