

Product Datasheet

Cancer Antigen 15-3 (CA 15-3) MUC 1 Antibody HM904

Product Name	Cancer Antigen 15-3 (CA 15-3) MUC 1 Antibody
Catalog Number	HM904
Description	Purified Monoclonal Antibody
Species/Host	CHO Cells
Tested Applications	ELISA, CLIA
Immunogen	Recombinant CA15-3
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffer, pH 7.4 with 0.15M NaCl
Specificity	CA 15-3 carbohydrate
Preservatives	0.1% Proclin 300