

Product Datasheet

CYFRA21-1 (Cytokeratin Fragment) Antibody HM908

Product Name	CYFRA21-1 (Cytokeratin Fragment) Antibody
Catalog Number	HM908
Description	Purified Monoclonal Antibody
Species/Host	CHO Cells
Tested Applications	ELISA, CLIA
Immunogen	Recombinant CYFRA21-1
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffer, pH 7.4 with 150 mM NaCl
Matched Pairs	HM907 (capture)
Preservatives	0.1% Proclin 300