

## Product Datasheet

### Inhibin B (INHB) Antibody HM910

<b>Product Name</b>	Inhibin B (INHB) Antibody
<b>Catalog Number</b>	HM910
<b>Description</b>	Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA
<b>Immunogen</b>	Recombinant INHB
<b>Alternative Names</b>	Inhibin-B
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	-20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 90 %
<b>Buffer</b>	10 mM Phosphate Buffer, pH 7.4 with 150 mM NaCl
<b>Specificity</b>	Inhibin B a subunit
<b>Matched Pairs</b>	HM909 (capture)
<b>Preservatives</b>	0.1% Proclin 300