

Product Datasheet

Human Immunoglobulin M (IgM) Antibody HM915

Product Name	Human Immunoglobulin M (IgM) Antibody
Catalog Number	HM915
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Human IgM
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with human IgA, IgD, IgE, IgG
Specificity	Human IgM
Preservatives	0.09% Sodium Azide