

Product Datasheet

Hepatitis B Surface Antigen (HBsAg) Antibody HM917

Product Name	Hepatitis B Surface Antigen (HBsAg) Antibody
Catalog Number	HM917
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	HBsAg
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Specificity	HBsAg, subtypes ad & ay
Matched Pairs	HM919 (detection)
Preservatives	0.09% Sodium Azide