

Product Datasheet

Hemoglobin (Hb) Antibody HM922

Product Name	Hemoglobin (Hb) Antibody
Catalog Number	HM922
Species/Host	Mouse
Tested Applications	LF, ELISA
Immunogen	Human hemoglobin
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with pig, bovine, horse and rabbit
Specificity	A0, A2, S and methemaglobin
Matched Pairs	HM923 (capture)
Preservatives	0.09% Sodium Azide