

Product Datasheet

Chlamydia pneumoniae Antibody HM925

Product Name	Chlamydia pneumoniae Antibody
Catalog Number	HM925
Species/Host	Mouse
Tested Applications	ELISA, IFA
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	No reactivity with C. trachomatis, K. pneumoniae, M. pneumoniae, L. pneumophila, S. pneumoniae, Influenza A/B, Parainfluenza, RSV
Specificity	C. pneumoniae
Preservatives	0.1% Sodium Azide