

Product Datasheet

Human IL-8 ELISA Set 851.530.015

Product Name	Human IL-8 ELISA Set
Catalog Number	851.530.015
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	31.25 -1000 pg/ml
Sensitivity	12.3 pg/ml
Tested Applications	ELISA
Alternative Names	CXCL8
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines
Specificity	Recognizes both natural and recombinant human IL-8
Sample Types	Serum, Plasma, Cell culture supernatant