

Product Datasheet

Phencyclidine (PCP) Antibody HM927

Product Name	Phencyclidine (PCP) Antibody
Catalog Number	HM927
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	PCP-KLH
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl
Cross Reactivity	Tenocyclidine 0.6%
Specificity	Phencyclidine
Matched Pairs	LA502 or LA018 (capture)
Preservatives	0.05% Sodium Azide