

Product Datasheet

Phencyclidine (PCP) Antibody HM929

| | |
|----------------------------|---------------------------------------------------------------------------|
| Product Name | Phencyclidine (PCP) Antibody |
| Catalog Number | HM929 |
| Species/Host | Mouse |
| Tested Applications | LF, ELISA |
| Immunogen | PCP-KLH |
| Brand | EastCoast Bio |
| Form/Appearance | Ion Exchange Purified Monoclonal Antibody |
| Storage | Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Purity | >90% by SDS-PAGE |
| Buffer | 0.015M Potassium Phosphate Buffer, pH 7.2 with 0.85% NaCl |
| Cross Reactivity | Tenocyclidine 10% |
| Specificity | Phencyclidine |
| Matched Pairs | LA018 (capture) |
| Preservatives | 0.05% Sodium Azide |