

Product Datasheet

Dengue Virus (DENV-NS1 Serotype 4) Antibody HM932

Product Name	Dengue Virus (DENV-NS1 Serotype 4) Antibody
Catalog Number	HM932
Description	Purified Monoclonal Antibody
Species/Host	CHO Cells
Tested Applications	LF
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 90 %
Buffer	10 mM Phosphate Buffered Saline, pH 7.4 with 150 mM NaCl
Specificity	NS1 protein of Dengue Virus type 4
Matched Pairs	HM513 (capture) & HM931 (detection)
Preservatives	0.1% Proclin 300