

## Product Datasheet

### Saint Louis Encephalitis Virus (SLEV-NS1) Antibody HM940

<b>Product Name</b>	Saint Louis Encephalitis Virus (SLEV-NS1) Antibody
<b>Catalog Number</b>	HM940
<b>Description</b>	Ion Exchange Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, LF
<b>Immunogen</b>	Recombinant SLEV NS1
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Ion Exchange Purified Monoclonal Antibody
<b>Storage</b>	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Isotype</b>	IgG2a
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 90 % by SDS-PAGE
<b>Buffer</b>	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
<b>Cross Reactivity</b>	No reactivity with Japanese Encephalitis Virus NS1 and Tick-Borne Encephalitis Virus NS1
<b>Specificity</b>	Saint Louis Encephalitis Virus
<b>Matched Pairs</b>	HM941, HM942 or HM943 (detection)
<b>Preservatives</b>	0.05% Sodium Azide