

Product Datasheet

Saint Louis Encephalitis Virus (SLEV-NS1) Antibody HM944

Product Name	Saint Louis Encephalitis Virus (SLEV-NS1) Antibody
Catalog Number	HM944
Description	Ion Exchange Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Recombinant SLEV NS1
Brand	EastCoast Bio
Form/Appearance	Ion Exchange Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG2b
Clonality	Monoclonal
Purity	> 90 % by SDS-PAGE
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Cross Reactivity	No reactivity with Japanese Encephalitis Virus NS1 and Tick-Borne Encephalitis Virus NS1
Specificity	Saint Louis Encephalitis Virus
Matched Pairs	HM941, HM942 or HM943 (detection)
Preservatives	0.05% Sodium Azide