

Product Datasheet

Respiratory Syncytial Virus (RSV) Antibody HM947

Product Name	Respiratory Syncytial Virus (RSV) Antibody
Catalog Number	HM947
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Purified native RSV
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with Influenza A, Influenza B, and Adenoviruses
Specificity	RSV fusion protein
Matched Pairs	HM947 (self-pairs)
Preservatives	0.09% Sodium Azide