

Product Datasheet

Adenovirus Antibody HM948

Product Name	Adenovirus Antibody
Catalog Number	HM948
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Purified Native Adenovirus
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	+2-8 °C
Isotype	IgG1
Clonality	Monoclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with Influenza A, Influenza B, RSV, Rotavirus, Astrovirus, Norovirus GI/GII, Enterovirus
Specificity	Adenovirus
Matched Pairs	HM948 (self-pairs)
Preservatives	0.09% Sodium Azide