

Product Datasheet

Legionella pneumophila Antibody HM949

Product Name	Legionella pneumophila Antibody
Catalog Number	HM949
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Purified Native L. pneumophila
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG3
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4 ± 0.2
Cross Reactivity	No reactivity with S. pneumoniae
Specificity	L. pneumophila
Matched Pairs	HM950 (detection)
Preservatives	0.09% Sodium Azide