

Product Datasheet

Legionella pneumophila Antibody HM950

| | |
|----------------------------|--|
| Product Name | Legionella pneumophila Antibody |
| Catalog Number | HM950 |
| Description | Purified Monoclonal Antibody |
| Species/Host | Mouse |
| Tested Applications | ELISA, LF |
| Brand | EastCoast Bio |
| Form/Appearance | Purified Monoclonal Antibody |
| Storage | Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing. |
| Isotype | IgG2b |
| Clonality | Monoclonal |
| Purity | ≥ 95 % by SDS-PAGE |
| Buffer | 10 mM Phosphate Buffered Saline, pH 7.4 |
| Cross Reactivity | No reactivity with S. pneumoniae |
| Specificity | L. pneumophila |
| Matched Pairs | HM949 (capture) |
| Preservatives | 0.09% Sodium Azide |