

## Product Datasheet

### Clostridium difficile toxin A Antibody HM951

<b>Product Name</b>	Clostridium difficile toxin A Antibody
<b>Catalog Number</b>	HM951
<b>Description</b>	Purified Monoclonal Antibody
<b>Species/Host</b>	Mouse
<b>Tested Applications</b>	ELISA, LF
<b>Immunogen</b>	Purified native C. difficile toxin A
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Monoclonal Antibody
<b>Storage</b>	+2-8 °C
<b>Isotype</b>	IgG1
<b>Clonality</b>	Monoclonal
<b>Purity</b>	> 95 % by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Cross Reactivity</b>	No reactivity with C. difficile toxin B, C. difficile GDH
<b>Specificity</b>	C. difficile toxin A
<b>Matched Pairs</b>	HM952 (capture)
<b>Preservatives</b>	0.09% Sodium Azide