

# Product Datasheet

## Human IL-10 ELISA Set 851.540.015

Product Name	Human IL-10 ELISA Set
Catalog Number	851.540.015
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	12.5-400 pg/ml
Sensitivity	5 pg/ml
Tested Applications	ELISA
Alternative Names	NULL
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human cytokines. Cross reactivity with simian IL-10. No cross reactivity with viral IL-10. No cross reactivity with viral IL-10.
Specificity	Recognizes both natural and recombinant human IL-10
Sample Types	Serum, Plasma, Cell culture supernatant