

Product Datasheet

Clostridium difficile toxin B Antibody HM954

Product Name	Clostridium difficile toxin B Antibody
Catalog Number	HM954
Description	Purified Monoclonal Antibody
Species/Host	Mouse
Tested Applications	ELISA, LF
Immunogen	Purified native C. difficile toxin B
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Isotype	IgG1
Clonality	Monoclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Cross Reactivity	No reactivity with C. difficile toxin A
Specificity	C. difficile toxin B
Matched Pairs	HM953 (capture)
Preservatives	0.09% Sodium Azide