

Product Datasheet

Borrelia burgdorferi Antibody J101

Product Name	Borrelia burgdorferi Antibody
Catalog Number	J101
Species/Host	Rabbit
Tested Applications	ELISA, WB, IHC
Immunogen	B. burgdorferi whole cells
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Cross reacts with Treponema pallidum, B. hermsii, and B. parkerii
Specificity	Reactive with WB bands present at 83kD, 41kD, 34kD, 31kD and other low molecular weight bands.
Preservatives	0.1% Sodium Azide