

Product Datasheet

Cytomegalovirus (CMV) Antibody J105

Product Name	Cytomegalovirus (CMV) Antibody
Catalog Number	J105
Species/Host	Goat
Tested Applications	ELISA, WB, IFA
Immunogen	Purified CMV strain AD169
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline pH 7.2
Cross Reactivity	Does not cross-react with HSV or VZV
Specificity	AD 169
Preservatives	0.1% Sodium Azide