

Product Datasheet

Feline Leukemia Virus (p27) Antibody J106

Product Name	Feline Leukemia Virus (p27) Antibody
Catalog Number	J106
Species/Host	Goat
Tested Applications	ELISA, WB, IFA
Immunogen	Purified FeLV p27
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffer, pH 7.2
Specificity	p27
Preservatives	0.1% Sodium Azide